
PLA BIACORE DATA ACQUISITION MODULE 1.1.0

Release Notes

June 6, 2016

COPYRIGHT

©2016 by Stegmann Systems GmbH, Rodgau, Germany.
All rights reserved.

CONTACT

Stegmann Systems GmbH
Raiffeisenstr. 2 // C1, C2
63110 Rodgau
Germany

Phone: +49 6106 77010-0
Fax: +49 6106 77010-190

www.stegmannsystems.com
support@stegmannsystems.com

TABLE OF CONTENTS

Copyright.....	2
Contact.....	2
Table of Contents.....	3
Welcome	4
Release Dates.....	5
Release Notes.....	6
Version 1.1.0.12	6
New	6
Changed	6
Fixed.....	6
Known Issues.....	6
Version 1.0.0.7	6
New	6
Changed	6
Fixed.....	6
Known Issues.....	6

WELCOME

Welcome to the Release Notes of **PLA Biacore Data Acquisition Module 1.1.0**. This document contains the cumulated Release Notes showing all new features, changes and fixes for all released versions.

This document is primarily intended for Functional and System Administrators.

RELEASE DATES

Product	Version	Release Date
PLA Biacore Data Acquisition Module	1.0.0.7	02/04/2016
	1.1.0.12	06/06/2016

Versions not listed here are internal versions and have never been released to the public. The release notes document the changes between the listed versions.

RELEASE NOTES

VERSION 1.1.0.12

Service release for PLA 3.0.4 compatibility.

NEW

- Compatibility with PLA 3.0.4.
- New import strategy "Start writing at first row (resets previously imported values)".

CHANGED

- N/A

FIXED

- Values were imported from wrong column.

KNOWN ISSUES

- N/A

VERSION 1.0.0.7

Initial release for PLA 3.0.3.

NEW

- Import of values from Biacore XML files.

CHANGED

- N/A

FIXED

- N/A

KNOWN ISSUES

- N/A