



MICROBIAL ASSAYS FOR ANTIBIOTICS

VERSION 1.0.0 (BUILD 45)

RELEASE NOTES

Version 1.0.0 Release Notes
REV 1.0

support@bioassay.de
+49 6106 77010-0
February 22, 2022

WELCOME

Welcome to the Release Notes of the Microbial Assays for Antibiotics Package Version 1.0.0 (Build 54). This document contains the cumulated Release Notes showing all new features, changes and fixes for all versions of the Microbial Assays for Antibiotics Package since its initial release.

This document is primarily intended for Functional and System Administrators.

TABLE OF CONTENTS

Welcome.....	2
Microbial Assays for Antibiotics Package Version 1.0.0 (Build 45)	4
New Features.....	4

MICROBIAL ASSAYS FOR ANTIBIOTICS PACKAGE VERSION 1.0.0 (BUILD 45)

New Features

- ◆ Initial release of the “Cylinder plate assay (USP <81>)” document type for performing cylinder plate assays in accordance with United States Pharmacopeia Chapter <81>.
- ◆ Initial release of the “Cylinder plate assay (JP 4.02)” document type for performing cylinder plate assays in accordance with Japanese Pharmacopoeia Chapter 4.02.
- ◆ Initial release of the “Turbidimetric assay (USP <81>/JP4.02)” document type for performing turbidimetric assays in accordance with United States Pharmacopeia Chapter <81> and Japanese Pharmacopoeia Chapter 4.02.
- ◆ Initial release of the “Combination of microbial assays” document type for performing combination of assays calculations and determining confidence limits in accordance with United States Pharmacopeia Chapter <81>.